



### Adpt-OM32-A32 Adaptor Map

NOTE: This diagram shows the channel mapping within the adaptor itself, and is not a map of the final pin-out.




### OM32 Connector



### BrainAmp Connector

UPPER ROW G 31 30 G 27 26 G 23 22 G 19 18 R 15 14 G 11 10 G 7 6 G 3 2 G



LOWER ROW G 32 G 29 28 G 25 24 G 21 20 G 17 16 G 13 12 G 9 8 G 5 4 G 1

G = Ground  
R = Reference